## How to Use Aerolizer

Patient Education Guide



Learn with - An Interactive training - Engaging animations - Your web cam Good Inhaler Technique -Increases medicine effect - Saves of cost of medicine - Reduces ER Visits

## www.use-inhalers.com

- For cleaning refer to the patient instructions that accompanies your inhaler.
- Following these step by step instructions will help you to inhale medicine effectively.



**Step1:** Remove the Aerolizer cover.



**Step2:** Hold the base of the Aerolizer firmly. Rotate the mouthpiece in the direction of the arrow to open it.



**Step3:** Insert the capsule in the base chamber of the Aerolizer.



**Step4:** Rotate the inhaler's mouthpiece to close it.



Step5: Press both the buttons once. You'll hear a click. Doing so, makes holes in the capsule and easily releases the medicine, when you are breathing in.



**Step6:** Keep your head upright. Hold the inhaler's mouthpiece horizontally. The buttons should be on the left and right side of the mouthpiece.



**Step7:** Exhale away from the inhaler.



**Step8:** Keep the mouthpiece in your mouth and take a quick and deep breath. When you do this, you'll hear the capsule vibrate.



**Step9:** Hold your breath for 10 seconds. Exhale slowly through your mouth or nose.



**Step10:** Take another breath in through the mouthpiece. This ensures that you get your complete dose of medicine.



**Step11:** Open the inhaler's mouthpiece, remove the used capsule, and throw it away. Do not keep capsules in the inhaler.



**Step12:** Close the mouthpiece and replace the cover.

Tip: Remember to rinse your mouth thoroughly after taking your dose and spit out the water.

Use-inhalers website is a great resource for Asthma and COPD patients. Learn how to use your Inhaler with interactive audio-visual instructions. Practice using the innovative webcam feature.

The content of this document is the property of use-inhalers.com. Use-inhalers.com owns all rights to the content. You agree not to copy any content or feature there in.